KULLANMA KILAVUZU

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INSTRUCTION FOR USE

1. Product Description

VIDA prostheses are used to form connections between abutments and copings. VIDA dental implant abutment systems consist of different types, sizes and platforms.

2. Intended use

VIDA dental implant abutment systems;

It is designed to be used as an implant abutment added to support artificial dental prostheses used to restore chewing function.

3. Indications

- Complete and partially edentulous patients
- In chin and facial defects
- Complete tooth deficiencies characterized by excessively resorbed ridge
- In patients who have difficulties in the use of removable partial dentures
- In patients who refuse to use removable prosthesis
- In patients who want to have their teeth prepared
- In fixed prosthesis patients with very long spaces
- It can be used even in individuals with excessive overbite, bruxism or foreign body biting habits, due to its high fracture strength, in any edentulous area or in severe damage to the soft tissues where the prosthesis sits.
- In cases where the height of the abutment is less than 7 mm and the axial thickness is less than 0.7 mm due to the patient's mouth closure, and in cases where the abutment must be angled more than 30 times surgically, it is more used than ceramic abutments in the posterior region.
- They can be used in the posterior shank and as a bridge abutment. Metal abutments themselves have no risk of breakage.

4. Contraindication

- Active Infection
- Foreign body sensitivity / material allergy
- The condition of the bone structure and insufficient bone density
- Unstable
- Relaxation
- Dislocation
- When the edentulous space is narrowed by the opposing and adjacent teeth
- In cases where the crown length of the abutment teeth is too short,
- In the presence of parafunctional habits such as bruxism
- When the use of a winged bridge (Kantilever) is designed
- Auxiliary teeth are not used if they lack adequate periodontal support.
- Your surgeon; failure to follow surgical operation prescription and/or inadequate surgical practice
- Inability of the edge of the cornea to adapt to the scallop structure of the gingiva;
- Deep gingival pockets formed as a result of not placing the implant deeper in order to provide a suitable aesthetic,
- Difficulty in cementation of the crown and removal of cement residues,
- The reflection in metallic blue color visible under the gingiva, Aesthetic negativity in patients with thin gingival structure or high smile lines,
- When the patient smiles, the arm of the abutment will be seen as a metal band, worse than the metal reflection from the arm. This metal tape poses an aesthetically unacceptable problem.

5. Warnings and Precautions

Before the operation, the patient should be informed about the surgical risks and the positive and negative effects should be explained.

-The patient should be warned that the implanted device has a certain lifespan during which normal bone cannot replace, may be damaged by strenuous activity or trauma, and may be replaced in the future.

- It is recommended to review the product-specific surgical technique before the surgeon performs the surgical procedure. VIDA DENTAL can provide surgical technical information. Please contact VIDA DENTAL sales representative.

- Choosing the right implant is extremely important. Appropriate type, size, largest loss component is taken into consideration, patient age, activity levels, bone density are preferred depending on whether or not he has had any previous surgical operation.

- VIDA DENTAL prostheses should be applied by securing respiration during oral use. Components or tools may injure the patient.

-Failure to follow these instructions may result in any or all of the following;

Adverse respiration

Ingestion of an ingredient

Prolonged healing time

- VIDA DENTAL abutments and copings have not been evaluated for safety and compatibility in a magnetic resonance environment.

- VIDA DENTAL abutments and copings have not been evaluated in terms of heating and collapse in a magnetic resonance environment.
- Small diameter abutments should not be used where high mechanical loads are required..
- Do not use implants made of different metals and supplied from different manufacturers together.
- The product is for single use. Do not use again.
- Check the expiration date of the product before use.
- Do not use the product if a foreign substance or contamination is seen in the package or on the product.
- Do not use the product that has fallen on the ground.
- Since there is a risk of infection, dispose of the product as medical waste after use within the framework of legal procedures. VIDA implant abutment systems and prosthetic parts are available in a variety of configurations. Abbreviations are used on each product label to help. Remove the healer cap or temporary healer and rinse thoroughly then dry the inside of the abutment. Remove the optimizer from the working model.
- After the abutment is placed properly, the screw is tightened. Be sure to fix the abutment with the appropriate screw.
- Torque Usage; When torques are greater than 35 Ncm, failure of the abutment and/or implants may occur. A torque value less than the recommended value may cause the abutment to loosen. It may lead to abutment and implant loss.

Selection of the product

In the selection of abutment superstructure, for the implant with a hexagonal implant connection of 2.40 and an internal connection angle of 22°, "2.48 Hexagon 22° B." Superstructures starting with

"2.40 Hexagon 89° B" for the implant with a hexagonal 2.40 connection and an 89° internal connection. Superstructures starting with, "3.10 Octagon B" for the implant whose implant connection is octagonal and 3.10. It is necessary to use abutments starting with "2.40 aesthetic" for superstructures and implants with a hexagonal 2.40 implant connection. However, in use, the abutment attached to the implant is checked and placed by the expert and when there is any rotation, the product is removed.

6. User Audience:

The person or persons who will use our device;

- Professionals who have completed all their medical education
- Have serious operational experience
- Be careful during the operation
- Must have sufficient experience in using the device

7. Patient Population:

Implants may not be preferred in young individuals who have not completed their development. In young individuals who have not reached the age of 17-18 in men and 16-17 in women, a decision is made whether to start the treatment by calculating the development and predicting the development before making a screw tooth. However, the development of each individual differs from each other. Therefore, the age ranges stated here are general averages.

- In patients who have not completed their Bone Development;

Since bone development continues in individuals who have not completed puberty, they are not suitable for implant application. In order to make an implant, bone development must be completed.

- In diabetic patients;

Implantation in uncontrolled diabetes patients is among the risky groups. Although tissue healing is slow in diabetic patients, the risk of infection increases.

- In patients with high cigarette consumption;

Smoking increases the risk of infection by preparing the ground for the bacterial plaque formed in the tissues in the mouth. The negative effect of cigarette on the integration of the implant with the bone due to carbon monoxide and toxic effects.

- Those with hypertension;

In patients with hypertension, there is an overreaction in situations such as pain and stress. In case of overreaction to stimuli, blood pressure may rise acutely during dental treatment and complications such as congestive heart failure or, rarely, stroke, bleeding may occur.

- Those with heart disease;

Implants should not be applied without consulting the doctor.

- Blood Coagulation Disorders;

Implants are not applied to those with blood clotting disease.

8. Shelf life:

The product does not have a specified shelf life. The lifespan of the products has been determined as 10 years based on literature support.

9. Compatibility Information

MR Compatibility Information

Our company has not conducted any tests regarding MRI Compatibility. MRI compatibility is unknown. Please inform your doctor or radiologist before having an MRI scan. Notes. The raw material of the product is Titanium.

10. Cleaning and Disinfecting

VIDA abutment and prosthesis parts are delivered in a non-sterile state. Products must be cleaned, disinfected and sterilized before being placed in the patient's mouth. Before using the devices, the recommendations in the procedure should be followed and cleaned, disinfected and sterilized;

- It should be washed under running water by brushing the inside and outside with a brush.
- You can clean the products manually using ultrasonic support, automatic method or disinfection methods as pre-treatment.
- When using the automatic method and disinfection method, the appropriate detergent should be selected and the manufacturers' instructions for use should be followed.

11. Sterilization

NOTE: Sterilization materials that are used by users you are sure of and have appropriate biological indicators should be used.

Our company's products are offered to the market as non-sterile, and the following method must be followed before sterilization.

Recommended Sterilization Method: The product should be removed from its current packaging and packaged with a packaging material suitable for autoclave sterilization.

System elements are delivered in clean condition, but they are not sterile. Before use, reusable hand tools must be pre-cleaned and all system elements must be sterilized by following the methods specified below.

In accordance with ISO 17665-1:2006, AAMI TIR 12:2004 and other relevant standards

METHOD	CYCLE	TEMPERATURE	MINIMUM DURATION
VAPOR	Vacuum	121 °C (273 °F)	Min. 30 minutes

TABLE-1

12. Instructions for Use

12.1 Procedure

Be sure to protect the prosthetic attachment by securing the abutment during polishing or other laboratory procedures.

- Insert the abutment into the implant analog on the study model.
- Make sure the analog-abutment connection holder elements are properly aligned.
- Straighten the abutment by tightening the main screw by hand.

12.2 Technical

The original abutment is taken from the laboratory. Remove the healing cap or temporary healer and rinse thoroughly, then dry the inside of the implant. Remove the optimizer from the working model.

- Cleaning, disinfection and sterilization of dental healers are described in Chapters 7 and 8.
- The healer is placed in the patients' mouth.

NOTE: After the abutment is placed properly, the screw is tightened.

- Make sure to fix the abutment implant with the appropriate screw.
- Tighten the abutment screw with the relevant screwdriver.

12.3 Torque Usage

When torques are greater than 35 Ncm, failure of the abutment and/or implants may occur. A torque value less than the recommended value may cause the abutment to loosen. It may lead to abutment and implant loss.

NOTE: The abutment or its components are tightened to the specified torque. The screws must not be removed, otherwise a new screw must be used.

Follow these instructions for the safe and correct use of VIDA products. It is the user's responsibility to use the device in accordance with the instructions.

The VIDA product is used as part or connection of an overall concept. Product sales made by third parties are not VIDA distributors. It does not contain any warranty or liability.

12.4 Usage areas of products

- Finger switch; Removes the healing cap inside the mouth.
- Gingiva Former; It remains in the patient's mouth for 1 week and 10 days, allowing the patient's gingiva to take shape.
- Measurement Piece; It ensures that the shape of the patient's mouth and the implant are exactly outside, allowing the technician to work on the model.
- Abutment; Implant abutments remain in the mouth for a lifetime after the healing period of the patient's implant in the mouth is completed, so that a crown can be placed on it that is both aesthetic and functional.
- Analog; It is the part used by the technician to measure the implant size in the mouth and work on the model.
- Abutment Screw; It is the part used to fix the implant and abutment.

13. Storage and Shipment Conditions

The products are stored at room temperature before being shipped to the customer. There are no special storage conditions.

14. Compatibility Table

Brand / model information of the implants that the dental abutment is compatible with are given in the table below.

Compatibility Table	
Code	Brand
BR	Bredent®
CM	Cowellmedi®
ZM	Zimmer®
EBL-ETS	EBI®
NB	NOBEL Biocare®
OSS-OTS	OSSTEM®
NEO	Neo Biotech®
DT	Dentium®
DS	Dentis®
BG	BEGO®
STS-SBL	Straumann®
ZD	Zinedent-Neodent®
MG	Megagen®

TABLE-2

15. Accessing the User Manual and Requesting a Hard Copy

The user manual is presented digitally in the form of a QR code on the product label. The user manual can be accessed by scanning the QR code with the camera applications of iOS and Android devices. In case the QR code is corrupted or cannot be read, the access link is located below the QR code.

If there is any problem in accessing the user manual over the internet and the situation is reported to Vida Dental, the printed version of the user manual will be delivered to the user within 7 days.

See e-IFU: www.vidadental.com.tr/ifu2.pdf

16. Symbols and Their Meanings Prepared in Accordance with EN 15223-1

	Notified Body Number		Non-Sterile
	Reference Number		Lot Number
	Do not use twice		Production Date
	Check Instructions for use		Must be used by experts
	Manufacturer Information		Barcode Number
	Warnings		Firm Logo
	Do not use if packaging is damaged		Keep away from direct sunlight
	Keep Dry		Medical Device
	UDI-DI		Do not sterilize twice
	Instructions for use QR code		

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